

```

Multi-transmit = "no";
Coil selection 1 = "SENSE-Head-32P";
    element selection = "selection 1";
Coil selection 2 = "SENSE-Head-32AH";
    element selection = "selection 1";
Dual coil = "yes";
CLEAR = "yes";
    body tuned = "no";
FOV      RL (mm) = 240;
        AP (mm) = 240;
        FH (mm) = 120;
Voxel size  RL (mm) = 3;
        AP (mm) = 3;
Slice thickness (mm) = 3;
Recon voxel size (mm) = 3;
RFOV (%) = 100;
Small FOV imaging = "no";
Fold-over suppression = "no";
Matrix scan = 80;
    reconstruction = 80;
Scan percentage (%) = 100;
SENSE = "yes";
    P reduction (AP) = 2;
    P os factor = 1;
k-t BLAST = "no";
Stacks = 1;
    type = "parallel";
    slices = 40;
    slice thickness (mm) = 3;
    slice gap = "user defined";
        gap (mm) = 0;
    slice orientation = "transverse";
    fold-over direction = "AP";
    fat shift direction = "P";
Stack Offc. AP (P=+mm) = 1.75740623;
        RL (L=+mm) = -0.521792829;
        FH (H=+mm) = -6.37945843;
    Ang. AP (deg) = 3.93005323;
        RL (deg) = -16.8998699;
        FH (deg) = 3.35096574;
Minimum number of packages = 1;
Slice scan order = "default";
Large table movement = "no";
PlanAlign = "no";
REST slabs = 0;
Interactive positioning = "no";
Patient position = "head first";
    orientation = "supine";
Scan type = "Imaging";
Scan mode = "MS";

```

```

    technique = "SE";
Modified SE = "no";
Acquisition mode = "cartesian";
Fast Imaging mode = "EPI";
    shot mode = "single-shot";
Echoes = 1;
    partial echo = "no";
TE = "user defined";
    (ms) = 80;
Flip angle (deg) = 90;
TR = "shortest";
Halfscan = "yes";
    factor = 0.600000024;
Water-fat shift = "minimum";
Shim = "auto";
Fat suppression = "SPIR";
    strength = "strong";
    frequency offset = "default";
Water suppression = "no";
BB pulse = "no";
MTC = "no";
Research prepulse = "no";
Diffusion mode = "DTI";
    sequence = "SE";
    gradient duration = "minimum";
    gradient overplus = "yes";
    directional resolution = "high";
    nr of b-factors = 4;
    b-factor order = "user defined";
    b-factors = 0, 1000, 1500,
    2000, (12) 0;
    average high b = "no";
SAR mode = "high";
B1 mode = "default";
PNS mode = "moderate";
Gradient mode = "default";
SoftTone mode = "no";
Cardiac synchronization = "no";
Respiratory compensation = "no";
Navigator respiratory comp = "no";
Flow compensation = "no";
Temporal slice spacing = "default";
NSA = 1;
Manual start = "no";
Dynamic study = "no";
    dyn stabilization = "yes";
Arterial Spin labeling = "no";
Preparation phases = "auto";
Interactive F0 = "no";
B1 field map = "no";

```

MIP/MPR = "no";  
Images = " M", (3) " no";  
Autoview image = " M";  
Calculated images = (4) " no";  
Reference tissue = "White matter";  
Preset window contrast = "soft";  
Reconstruction mode = "immediate";  
Save raw data = "no";  
Hardcopy protocol = "no";  
Ringing filtering = "default";  
Geometry correction = "default";  
IF\_info\_seperator = 1634755923;  
Total scan duration = "10:48.1";  
Rel. signal level (%) = 100;  
Act. TR (ms) = "4577";  
Act. TE (ms) = "80";  
ACQ matrix M x P = "80 x 78";  
ACQ voxel MPS (mm) = "3.00 / 3.07 / 3.00";  
REC voxel MPS (mm) = "3.00 / 3.00 / 3.00";  
Scan percentage (%) = 97.727272;  
Packages = 1;  
Min. slice gap (mm) = 0;  
Diffusion gradient timing DELTA / delta (ms) =  
"51.1 / 18.0";  
EPI factor = 43;  
WFS (pix) / BW (Hz) = "11.334 / 38.3";  
BW in EPI freq. dir. (Hz) = "2668.1";  
SAR / head = "< 33 %";  
Whole body / level = "< 0.1 W/kg / normal";  
B1 rms = "1.33 uT";  
PNS / level = "65 % / normal";  
Sound Pressure Level (dB) = 15.4141731;