

TRACOE VARIO **PRODUCT LINE**

TRACOE® vario

TRACOE vario product line

In most of the regions of the world double-lumen tubes like the TRACOE *twist* line are the favourable choice of tube. But in some areas the TRACOE *vario* line with single lumen has become more and more popular due to the need for XL products.

TRACOE offers 3 different types of *vario* tubes (spiral-reinforced, transparent and with subglottic suction line) available in standard and extra-long.

The soft and flexible single lumen TRACOE *vario* tracheostomy tubes have a patented variably adjustable neck flange. This allows for a customisable adjustment of the tube suitable for different anatomies and people with special needs. The flange can be moved along the tube by using the simple and easy self-securing push-button mechanism. For additional security, the neck flange can be securely locked by shifting up the orange lever system. The flexible wings on the flange can be rotated independently. This gives maximum flexibility to accommodate even the most delicate and difficult anatomical conditions in respect of the tracheostoma.

Features	Benefits
Patented adjustable neck flange	perfect adjustment to the anatomy of the patient's neck
	 easy to handle self-securing push-button
	no screw mechanism
	 safe and reliable locking system
Pliable wings	maximum flexibility to accommodate even the most delicate and difficult anatomies
	 individually adjustable and soft wings
	 neck straps can easily be attached while the tube is in place
Soft PVC cuff	Please note: Do not use cremes which contain panthenol (Dexpanthenol, Bepanthen! (hardens the material)

Features	Benefits
Spiral-reinforced tubes	 Due to its spiral reinforcement, the tube is kink-free and has great flexibility at the same time. The tube will not be pressed together or collapses. It is embedded within the soft tube material.
Transparent tubes	 Metal-free (suitable for MRI) Soft and shape retaining even under pressure (kink resistant) X-ray contrast line for easy and safe monitoring of the tube position in situ

A verification regarding the MRI ability of our tubes is separately available.

Type of Tube (for details please see Full Line Catalogue)	Clinical Function
REF 450 Tracheostomy Tube with Cuff, spiral reinforced	 Ventilation Protection from aspiration Provides safe airway until patient is weaned off of the ventilator
REF 451 Extra-long Tracheostomy Tube with Cuff, spiral reinforced	 Can be used for patients with strong necks and obese patients Extra length to bridge tracheal stenosis and/or injuries Ventilation / Protection from aspiration Provides safe airway until patient is weaned off of the ventilator

TRACOE vario

Type of Tube (for details please see Full Line Catalogue)	Clinical Function
REF 455 Tracheostomy Tube, spiral reinforced	 Patient does not need ventilation Patient can control his own secretions No significant risk of aspiration Access for suctioning secretions For stable patients during weaning process To bridge upper airway obstruction
REF 460 Transparent Tracheostomy Tube with Cuff	 Metal-free (suitable for MRI) Ventilation Protection from aspiration Provides safe airway until patient is weaned off of the ventilator
REF 461 Extra-long transparent Tracheostomy Tube with Cuff	 Can be used for patients with strong necks and obese patients Extra length to bridge tracheal stenosis and/or injuries Metal-free (suitable for MRI) Ventilation / Protection from aspiration Provides safe airway until patient is weaned off of the ventilator

Type of Tube (for details please see Full Line Catalogue)	Clinical Function
REF 470 Transparent Tracheostomy Tube with Cuff and integrated Subglottic Suction Line	Suction line is completely embedded into the tube wall
	 Reduced risk of aspiration and associated infection risk
	 Dysphagia management (swallowing difficulties)
USP = Only PVC tube on the market with integrated suction line and	 Suction port above the cuff allows for regular removal of secretions
adjustable neck flange!	Metal-free (suitable for MRI)
	 Ventilation
	 Provides safe airway until patient is weaned off of the ventilator
Please refer to the chapter on "Subglottic Suctioning"	Please note: In case the diameter of the tube is too large, it can happen that the suction port is blocked by the tracheal mucosa.
	 Can be used with strong necks and obese patients
REF 471 Extra-long transparent Tracheostomy Tube with Cuff and integrated Subglottic Suction Line	 Extra length to bridge tracheal stenosis and/or injuries
	 Suction line is completely embedded into the tube wall
	 Reduced risk of aspiration and associated infection risk
USP = Only PVC tube on the market with integrated suction line and adjustable neck flange!	 Dysphagia management (swallowing difficulties)
	 Suction port above the cuff allows for regular removal of secretions
	Metal-free (suitable for MRI)
	 Ventilation
Please refer to the chapter on "Subglottic Suctioning"	 Provides safe airway until patient is weaned off of the ventilator
	Please note: In case the diameter of the tube is too large, it can happen that the suction port is blocked by the tracheal mucosa.