

## >> Follows the movements of the tympanic membrane

- >> Moveable head plate counteracts tension
- >> Clinically evaluated CliP

## CliP PARTIAL FlexiBAL PROSTHESIS







CliP Partial FlexiBAL Material: Pure Titanium (ASTM F67 Medical Grade) Shaft Diameter: 0.2 mm

LENGTH L (mm)	FL (mm)	REF
1.75	0.75	1002 350
2.00	1.00	1002 351
2.25	1.25	1002 352
2.50	1.50	1002 353
2.75	1.75	1002 354
3.00	2.00	1002 355
3.25	2.25	1002 356
3.50	2.50	1002 357

FL: Functional Length

MRI information: www.kurzmed.com

Developed in close collaboration with the ENT clinics of the university hospitals in Cologne and Dresden, Germany.

Heinz Kurz GmbH Medizintechnik | Tuebinger Strasse 3 72144 Dusslingen | Germany | Phone: +49 (0)7072/9179-0 E-Mail: info@kurzmed.de | www.kurzmed.com

KURZ Medical, Inc. | 5126 South Royal Atlanta Drive Tucker, GA 30084 | USA Toll Free Order #: 866.449.8020 | Phone: 770.349.6330 | E-Mail: info@kurzmed.com | www.kurzmed.com

**( 6** 0124

Prescription Devices

Caution: United States Federal Law restricts devices identified herein to sale by or on the order of a physician. This brochure does not replace nor does it set forth the complete contents of the "Instructions for Use" for the products in this brochure, and is not a substitute for reviewing and understanding that important information. Therefore, before using any products included in this brochure, please review the entire contents of the respective "Instructions for Use" information.

Because of registration the products are currently not available in all countries

